


Anti-Arf1 Rabbit Polyclonal Antibody For Wb, Elisa - P84077

Item Number: CM0024098



Introduction

High-quality ARF1 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat ARF1. Suitable for Golgi trafficking, synaptic plasticity, and cholera toxin activation studies.

[Learn More](#)

Parameter	Value
Product Name	ARF1 Rabbit pAb
Remarks/Alias	PVNH8; ARF1
Species	Human
GeneID (Human)	375
GeneID	375
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human ARF1 (NP_001649.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84077
Protein Weight	21kDa
Shipping	Ice bag