

Anti-Arf12 Rabbit Polyclonal Antibody (Oryza Sativa Arf12) - Arf12

Item Number: CM0018999

Introduction

Anti-ARF12 Rabbit Polyclonal Antibody detects Oryza sativa ARF12 (Q0J951, 91kDa). Suitable for WB and ELISA. Validated in rice, this pAb is ideal for plant research and auxin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	ARF12 Rabbit pAb
Remarks/Alias	ARF8; ARF12; OsARF8; OsJ_16583
Species	Oryza sativa
Gene ID	4337363
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 710-818 of Oryza sativa ARF12. (Q0J951)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
Protein Weight	91kDa
Shipping	Ice bag