


Anti-Arf4 Polyclonal Antibody For Wb, If/Icc, Elisa - P18085

Item Number: CM0028620



Introduction

Rabbit polyclonal antibody against human ARF4 (ADP-ribosylation factor 4) validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Research use in Golgi trafficking and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	ARF4 Rabbit pAb
Remarks/Alias	ARF2; ARF4
Species	Human
Gene ID (Human)	378
Gene ID	378
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ARF4 (NP_001651.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18085
Protein Weight	21kDa
Shipping	Ice bag