

Anti-Arf5 (Ko Validated) Rabbit Polyclonal Antibody For Wb, Elisa - P84085

Item Number: CM0024118



Introduction

Rabbit polyclonal antibody against ARF5, validated for WB and ELISA in human, mouse, and rat. Ideal for studying ADP-ribosylation factor 5 in Golgi trafficking, mTORC1 signaling, and lysosomal quality control.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ARF5 Rabbit pAb
Remarks/Alias	ARF5; F5
Species	Human
GeneID (human)	381
GeneID	381
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ARF5 (NP_001653.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84085
Protein Weight	21kDa
Shipping	Ice bag