

Anti-Arf5 Monoclonal Antibody For Wb, If/Icc, Elisa - P84085

Item Number: CM0010122



Introduction

Rabbit monoclonal antibody against human ARF5, validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for studying Golgi trafficking and mTOR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ARF5 Rabbit mAb
Remarks/Alias	ARF5; ADP-ribosylation factor 5
Species	Human
Gene ID (Human)	381
Gene ID	381
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 81-180 of human ARF5 (P84085).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84085
Protein Weight	21 kDa
Shipping	Ice bag