

Anti-Arf5 Rabbit Polyclonal Antibody For Wb, Elisa - P84085

Item Number: CM0028662



Introduction

Rabbit polyclonal antibody targeting human ARF5. Validated for Western blot (WB) and ELISA. Detects endogenous ARF5 in human, mouse, and rat samples. Suitable for research on Golgi trafficking, mTORC1 signaling, and lysosomal quality control.

[Learn More](#)

Parameter	Value
Product Name	ARF5 Rabbit pAb
Remarks/Alias	ARF5
Species	Human
GeneID	381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ARF5 (NP_001653.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84085
Protein Weight	21kDa
Shipping	Ice bag