

## Anti-Arf6 Monoclonal Antibody For Wb, Elisa - P62330

Item Number: CM0003153



### Introduction

Rabbit monoclonal anti-ARF6 antibody validated for Western blot and ELISA, cross-reacting with human, mouse, and rat. Ideal for studying GTPase ARF6 in endocytic recycling, cytoskeletal dynamics, and neurodevelopment.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ARF6 Rabbit mAb                                                                                                     |
| Remarks/Alias    | ARF6; ADP-ribosylation factor 6                                                                                     |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 382                                                                                                                 |
| GeneID           | 382                                                                                                                 |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 76-175 of human ARF6 (P62330). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P62330                                                                                                              |
| Protein Weight   | 20kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |