

Anti-Arf6 Polyclonal Antibody For Wb, If/Icc, Elisa - P62330

Item Number: CM0014576



Introduction

Rabbit polyclonal antibody targeting human ARF6 (P62330), validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat ARF6; suitable for studying membrane trafficking and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	ARF6 Rabbit pAb
Remarks/Alias	ARF6
Species	Human
Gene ID (Human)	382
Gene ID	382
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human ARF6 (NP_001654.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62330
Protein Weight	20kDa
Shipping	Ice bag