


Anti-Arfgap1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N6T3

Item Number: CM0026369



Introduction

High-quality ARFGAP1 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat ARFGAP1. Ideal for membrane trafficking and Golgi dynamics research.

[Learn More](#)

Parameter	Value
Product Name	ARFGAP1 Rabbit pAb
Remarks/Alias	ARF1GAP; HRIHFB2281; ARFGAP1
Species	Human
GeneID (Human)	55738
GeneID	55738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ARFGAP1 (NP_001268411.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6T3
Protein Weight	45kDa
Shipping	Ice bag