


Anti-Arfgap2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N6H7

Item Number: CM0023487



Introduction

Rabbit polyclonal antibody against ARFGAP2, a GTPase-activating protein for ARF1 implicated in ER-Golgi transport. Validated for WB, IF/ICC, ELISA in human, mouse, rat samples.

[Learn More](#)

Parameter	Value
Product Name	ARFGAP2 Rabbit pAb
Remarks/Alias	IRZ; ZFP289; ZNF289; NBLA10535; ARFGAP2
Species	Human
GeneID (Human)	84364
GeneID	84364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 242-521 of human ARFGAP2 (NP_115765.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6H7
Protein Weight	57kDa
Shipping	Ice bag