


Anti-Arfgef1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6D6

Item Number: CM0019333



Introduction

Anti-ARFGEF1 Rabbit Polyclonal Antibody detects human ARFGEF1 (BIG1), a 209 kDa Golgi-associated guanine nucleotide-exchange factor. Suitable for WB and ELISA applications in vesicular trafficking research, ideal for studying Golgi structure, cell polarity, and integrin maturation in human samples.

[Learn More](#)

Parameter	Value
Product Name	ARFGEF1 Rabbit pAb
Remarks/Alias	BIG1; P200; DEDISB; ARFGEP1; ARFGEF1
Species	Human
GeneID (Human)	10565
GeneID	10565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human ARFGEF1 (NP_006412.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6D6
Protein Weight	209kDa
Shipping	Ice bag