

Anti-Arfrp1 / Arp Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13795

Item Number: CM0002588

Introduction

Rabbit monoclonal antibody targeting human ARFRP1/ARP (Q13795) for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying Golgi trafficking and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ARFRP1 / ARP Rabbit mAb
Remarks/Alias	ARP; Arp1; ARL18; ARFRP1 / ARP
Species	Human
GeneID (Human)	10139
GeneID	10139
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 102-201 of human ARFRP1 / ARP (Q13795)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13795
Protein Weight	23kDa
Shipping	Ice bag