

Anti-Arfrp1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13795

Item Number: CM0017599

Introduction

Anti-ARFRP1 rabbit polyclonal antibody for Western blot and ELISA. Recognizes human, mouse, and rat ARFRP1, a trans-Golgi GTPase regulating protein sorting and lipoprotein lipidation.

[Learn More](#)

Parameter	Value
Product Name	ARFRP1 Rabbit pAb
Remarks/Alias	ARP; Arp1; ARL18; ARFRP1
Species	Human
GeneID (Human)	10139
GeneID	10139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human ARFRP1 (NP_001254476.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13795
Protein Weight	23kDa
Shipping	Ice bag