


Anti-Arginase 1 (Arg1) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P05089

Item Number: CM0019168



Introduction

Rabbit polyclonal antibody targeting human Arginase 1 (ARG1), validated for Western blot, immunofluorescence/immunoperoxidase, and ELISA. Cross-reacts with mouse and rat. Suitable for studying urea cycle and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Arginase 1 (ARG1) Rabbit pAb
Remarks/Alias	ARG1; arginase-1; Arginase 1 (ARG1)
Species	Human
GeneID (human)	383
GeneID	383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human Arginase 1 (ARG1) (NP_000036.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05089
Protein Weight	35kDa
Shipping	Ice bag