

Anti-Arginase-1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q61176

Item Number: CM0019598

Introduction

Rabbit polyclonal antibody targeting mouse Arginase-1 (Arg-1), validated for WB, IF-P, and ELISA. Cross-reacts with rat. Ideal for urea cycle, macrophage polarization, and tumor immune evasion studies.

[Learn More](#)

Parameter	Value
Product Name	Arginase-1 Rabbit pAb
Remarks/Alias	AI; PGIF; Arg-1; Arginase-1
Species	Mouse
GeneID (Human)	383
GeneID (Mouse)	11846
Immunogen	Recombinant protein of mouse Arginase-1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q61176
Protein Weight	34KDa
Shipping	Ice bag