


Anti-Arginase 2 (Arg2) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P78540

Item Number: CM0003955



Introduction

Rabbit monoclonal antibody against human Arginase 2 (ARG2), validated for WB, IHC-P, ELISA. Targets mitochondrial arginase-2, a key regulator of arginine metabolism, immune response, and endothelial function. Ideal for immunology and metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	Arginase 2 (ARG2) Rabbit mAb
Remarks/Alias	ARG2; arginase 2; Arginase 2 (ARG2)
Species	Human
GeneID (human)	384
GeneID	384
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 255-354 of human Arginase 2 (ARG2) (P78540)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P78540
Protein Weight	39kDa
Shipping	Ice bag