

Anti-Argonaute-2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ukv8

Item Number: CM0004710



Introduction

KD Validated Argonaute-2 rabbit monoclonal antibody for Western blot, IHC-P, IF/ICC, ELISA. Validated in human, mouse, rat. Ideal for RNA interference and gene silencing studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Argonaute-2 Rabbit mAb
Remarks/Alias	PPD; Q10; CASC7; EIF2C2; LESKRES; LINC00980
Species	Human
GeneID (Human)	27161
GeneID	27161
Immunogen	A synthetic peptide corresponding to amino acids 101-200 of human Argonaute-2 (NP_036286.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKV8
Protein Weight	97kDa
Shipping	Ice bag