

Anti-Argonaute-2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ukv8

Item Number: CM0030556



Introduction

The Argonaute-2 Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, ELISA in human, mouse, rat. Argonaute-2 is a core RISC component mediating miRNA/siRNA-guided gene silencing for research.

[Learn More](#)

Parameter	Value
Product Name	Argonaute-2 Rabbit pAb
Remarks/Alias	PPD; Q10; CASC7; EIF2C2; LESKRES; LINC00980; Argonaute-2
Species	Human
GeneID (human)	27161
GeneID	27161
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Argonaute-2 (NP_036286.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKV8
Protein Weight	97kDa
Shipping	Ice bag