


Anti-Arhhgap1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q07960

Item Number: CM0027541



Introduction

High-quality rabbit polyclonal antibody against human ARHGAP1, a Rho GTPase-activating protein that regulates Cdc42, Rac, and Rho. Validated for WB, IP, and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP1 Rabbit pAb
Remarks/Alias	RHOGAP; RHOGAP1; CDC42GAP; p50rhoGAP; ARHGAP1
Species	Human
Gene ID (Human)	392
Gene ID	392
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-439 of human ARHGAP1 (NP_004299.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07960
Protein Weight	50kDa
Shipping	Ice bag