

Anti-Arhgap11B Rabbit Polyclonal Antibody For Wb, Elisa - Q3Krb8

Item Number: CM0020840

Introduction

High-quality ARHGAP11B Rabbit Polyclonal Antibody for Western blot and ELISA. Targets human ARHGAP11B (SWISS: Q3KRB8), a hominin-specific protein crucial for neocortex development. Ideal for neurological research.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP11B Rabbit pAb
Remarks/Alias	B ⁺ -T; FAM7B1; GAP (1-8); ARHGAP11B
Species	Human
GeneID (Human)	89839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-267 of human ARHGAP11B (NP_001034930.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q3KRB8
Protein Weight	30kDa
Shipping	Ice bag