

Anti-Arhgap22 Polyclonal Antibody For Wb, Elisa - Q7Z5H3

Item Number: CM0022001

Introduction

High-quality anti-ARHGAP22 polyclonal antibody validated for WB and ELISA. Detects human ARHGAP22, cross-reacts with mouse and rat. Ideal for angiogenesis and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP22 Rabbit pAb
Remarks/Alias	RhoGAP2; RhoGap22; ARHGAP22
Species	Human
GeneID (Human)	58504
GeneID	58504
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 510-650 of human ARHGAP22 (NP_001242953.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7Z5H3
Protein Weight	77kDa
Shipping	Ice bag