

Anti-Arhgap25 Rabbit Polyclonal Antibody For Wb, Elisa - P42331

Item Number: CM0017417



Introduction

Rabbit polyclonal antibody against ARHGAP25, validated for Western blot and ELISA. Detects human, mouse, and rat ARHGAP25. Suitable for Rho GTPase signaling studies. - P42331

[Learn More](#)

Parameter	Value
Product Name	ARHGAP25 Rabbit pAb
Remarks/Alias	KAIA0053; HEL-S-308; ARHGAP25
Species	Human
Gene ID (Human)	9938
Gene ID	9938
Immunogen	Recombinant fusion protein containing amino acids 1-370 of human ARHGAP25 (NP_001159748.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42331
Protein Weight	73kDa
Shipping	Ice bag