

Anti-Arhgap25 Rabbit Polyclonal Antibody For Wb, Ihc-P - P42331

Item Number: CM0017516



Introduction

Rabbit polyclonal antibody targeting human ARHGAP25 (Rho GTPase-activating protein 25), suitable for WB, IHC-P, IF/ICC, ELISA. Validated in human and mouse samples. Immunogen corresponds to amino acids 1-370.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP25 Rabbit pAb
Remarks/Alias	KAIA0053; HEL-S-308; ARHGAP25
Species	Human
Gene ID (Human)	9938
Gene ID	9938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-370 of human ARHGAP25 (NP_001159748.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42331
Protein Weight	73kDa
Shipping	Ice bag