

Anti-Arhgap4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P98171

Item Number: CM0022493



Introduction

Rabbit polyclonal antibody to ARHGAP4, a Rho GTPase-activating protein involved in stress fiber regulation. Validated for WB, IHC-P, IF/ICC, ELISA, reacts with human, mouse, rat species. Ideal for hematopoietic cell research.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP4 Rabbit pAb
Remarks/Alias	C1; RGC1; p115; SrGAP4; RhoGAP4; ARHGAP4
Species	Human
GeneID (Human)	393
GeneID	393
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 900 to the C-terminus of human ARHGAP4 (NP_001158213.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98171
Protein Weight	105kDa
Shipping	Ice bag