

Anti-Arhgap44 Polyclonal Antibody For Wb And Elisa - Q17R89

Item Number: CM0028291

Introduction

Rabbit polyclonal antibody targeting human ARHGAP44/RICH2 (Gene ID: 9912). Validated in WB and ELISA; cross-reacts with mouse and rat. Ideal for exploring Rho GTPase signaling, synaptic plasticity, and dendritic spine dynamics.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP44 Rabbit pAb
Remarks/Alias	RICH2; NPC-A-10; ARHGAP44
Species	Human
Gene ID (Human)	9912
Gene ID	9912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 666-772 of human ARHGAP44 (NP_055674.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q17R89
Protein Weight	89kDa
Shipping	Ice bag