

Anti-Arhgef10L Polyclonal Antibody (Rabbit) For Wb, Elisa - Q9Hce6

Item Number: CM0022695

Introduction

Rabbit polyclonal antibody against ARHGEF10L (GrinchGEF), validated for Western blot and ELISA in human and rat samples. Targets a Rho guanine nucleotide exchange factor involved in cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF10L Rabbit pAb
Remarks/Alias	GrinchGEF; ARHGEF10L
Species	Human
Gene ID (human)	55160
Gene ID	55160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1120-1230 of human ARHGEF10L (NP_001011722)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HCE6
Protein Weight	140kDa
Shipping	Ice bag