

Anti-Arhgef11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O15085

Item Number: CM0025909



Introduction

Anti-ARHGEF11 Rabbit Polyclonal Antibody suitable for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Detects human and mouse ARHGEF11, a Rho guanine nucleotide exchange factor involved in RhoA GTPase regulation and neurite outgrowth.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF11 Rabbit pAb
Remarks/Alias	GTRAP48; PDZ-RHOGEF; ARHGEF11
Species	Human
Gene ID (Human)	9826
Gene ID	9826
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ARHGEF11 (O15085).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15085
Protein Weight	168kDa
Shipping	Ice bag