

Anti-Arhgef3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr81

Item Number: CM0031007

Introduction

Rabbit polyclonal antibody targeting human ARHGEF3 (GEF3/STA3/XPLN), validated for WB and ELISA, detecting human and mouse species. Suitable for research on RhoA/RhoB signaling and cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF3 Rabbit pAb
Remarks/Alias	GEF3; STA3; XPLN; ARHGEF3
Species	Human
GeneID (Human)	50650
GeneID	50650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 287-526 of human ARHGEF3 (NP_062455.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NR81
Protein Weight	60kDa
Shipping	Ice bag