

Anti-Arhgef5 Polyclonal Antibody For Ihc-P, Elisa - Q12774

Item Number: CM0019465

Introduction

Rabbit polyclonal antibody against human ARHGEF5 (Ephexin-3), validated for IHC-P and ELISA. Cross-reacts with mouse and rat. Detects a 177 kDa protein. Suitable for studies of Rho GTPase signaling, actin dynamics, and oncogenic processes.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF5 Rabbit pAb
Remarks/Alias	P60, TIM, GEF5, TIM1, ARHGEF5
Species	Human
GeneID (human)	7984
GeneID	7984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1248-1597 of human ARHGEF5 (NP_005426.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12774
Protein Weight	177kDa
Shipping	Ice bag