

Anti-Arhgef7 Rabbit Polyclonal Antibody For Wb, Elisa - Q14155

Item Number: CM0013218



Introduction

Rabbit polyclonal antibody against human ARHGEF7 (Beta-Pix/COOL-1), validated for WB and ELISA. Cross-reacts with human and mouse; detects ~90 kDa Rho GEF protein involved in cell migration and synapse formation.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF7 Rabbit pAb
Remarks/Alias	P50; P85; PAK3; PIXB; COOL1; P50BP; COOL-1; P85SPR; BETA-PIX; P85COOL1; Nbla10314; ARHGEF7
Species	Human
Gene ID (Human)	8874
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ARHGEF7 (NP_001106984.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14155
Protein Weight	90kDa
Shipping	Ice bag