

Anti-Arhgef9 Rabbit Polyclonal Antibody For If/Icc And Elisa - O43307

Item Number: CM0021741



Introduction

Rabbit polyclonal antibody targeting human ARHGEF9 (Collybistin), validated for IF/ICC and ELISA. Recognizes endogenous ARHGEF9 protein with a predicted molecular weight of 61 kDa. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF9 Rabbit pAb
Remarks/Alias	DEE8; PEM2; EIEE8; PEM-2; HPEM-2; COLLYBISTIN
Species	Human
GeneID (Human)	23229
GeneID	23229
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43307
Protein Weight	61kDa
Shipping	Ice bag