

Anti-Arhgef9 Rabbit Polyclonal Antibody For Wb, Elisa - O43307

Item Number: CM0026943



Introduction

Rabbit polyclonal antibody against human ARHGEF9 (collybistin). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neuroscience and synaptic research.

[Learn More](#)

Parameter	Value
Product Name	ARHGEF9 Rabbit pAb
Remarks/Alias	DEE8; PEM2; EIEE8; PEM-2; HPEM-2; COLLYBISTIN; ARHGEF9
Species	Human
Gene ID	23229
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-516 of human ARHGEF9 (NP_056000.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43307
Protein Weight	61kDa
Shipping	Ice bag