

Anti-Arid1A Monoclonal Antibody For Ihc-P, If/Icc, Elisa - O14497

Item Number: CM0002468



Introduction

Rabbit monoclonal antibody against human ARID1A (SWISS: O14497). Suitable for IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Immunogen: recombinant protein fragment (aa 1204-1314). Ideal for chromatin remodeling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	ARID1A Rabbit mAb
Remarks/Alias	ELD; B120; CSS2; OSA1; P270; hELD; BM029; MRD14; hOSA1; BAF250; C1orf4; BAF250a; SMARCF1; ARID1A
Species	Human
Gene ID	8289 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1204-1314 of human ARID1A (O14497)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14497
Protein Weight	242kDa
Shipping	Ice bag