

Anti-Arid1A Rabbit Polyclonal Antibody For Wb, Elisa - O14497

Item Number: CM0023998



Introduction

ARID1A Rabbit pAb validated for WB and ELISA, reacting with human, mouse, and rat. Ideal for chromatin remodeling research in neural development and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ARID1A Rabbit pAb
Remarks/Alias	ELD; B120; CSS2; OSA1; P270; hELD; BM029; MRD14; hOSA1; BAF250; C1orf4; BAF250a; SMARCF1
Species	Human
GeneID (Human)	8289
GeneID	8289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2195-2285 of human ARID1A (NP_006006.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14497
Protein Weight	242kDa
Shipping	Ice bag