


Anti-Arid1B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nfd5

Item Number: CM0016408



Introduction

Rabbit polyclonal antibody against human ARID1B (BAF250B). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for chromatin remodeling and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ARID1B Rabbit pAb
Remarks/Alias	CSS1; OSA2; 6A3-5; DAN15; MRD12; P250R; BRIGHT; BAF250B; SMARCF2; ELD/OSA1; ARID1B
Species	Human
GeneID (Human)	57492
GeneID	57492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-650 of human ARID1B (NP_059989.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFD5
Protein Weight	244kDa
Shipping	Ice bag