


Anti-Arid2 Polyclonal Antibody For Wb, Ip, Elisa - Q68Cp9

Item Number: CM0025531



Introduction

Rabbit polyclonal antibody targeting human ARID2 (BAF200), a SWI/SNF chromatin remodeling complex subunit involved in transcriptional regulation. Validated for WB, IP, and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	ARID2 Rabbit pAb
Remarks/Alias	CSS6; p200; BAF200; SMARCF3; ARID2
Species	Human
GeneID (Human)	196528
GeneID	196528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1586-1835 of human ARID2 (NP_689854.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q68CP9
Protein Weight	197kDa
Shipping	Ice bag