

Anti-Arid5A Mouse Monoclonal Antibody For Wb, Elisa - Q03989

Item Number: CM0007679

Introduction

Anti-ARID5A mouse monoclonal antibody validated for WB and ELISA. Reacts with human and mouse ARID5A (64kDa). Ideal for studying transcriptional regulation, chondrocyte differentiation, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	ARID5a Mouse mAb
Remarks/Alias	MRF1; MRF-1; RP11-363D14
Species	Human
GeneID (Human)	10865
GeneID	10865
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-594 of human ARID5a (Q03989).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q03989
Protein Weight	64kDa
Shipping	Ice bag