

Anti-Arih1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y4X5

Item Number: CM0020845



Introduction

Rabbit polyclonal antibody against ARIH1, an E3 ubiquitin-protein ligase involved in ubiquitination and DNA damage response. Validated for WB and ELISA with human, mouse, and rat samples. Ideal for ubiquitin pathway studies.

[Learn More](#)

Parameter	Value
Product Name	ARIH1 Rabbit pAb
Remarks/Alias	ARI; HARI; HHARI; UBCH7BP; ARIH1
Species	Human
GeneID (Human)	25820
GeneID	25820
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400 to the C-terminus of human ARIH1 (NP_005735.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4X5
Protein Weight	64kDa
Shipping	Ice bag