


## Anti-Arl1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P40616

Item Number: CM0014505



### Introduction

Anti-ARL1 rabbit polyclonal antibody validated for western blot, immunofluorescence/ICC, and ELISA. Detects human, mouse, and rat ARL1, a Golgi-associated GTPase involved in vesicular trafficking. Ideal for basic research and diagnostic development.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | ARL1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | ARFL1; ARL1                                                                                                      |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 400                                                                                                              |
| Gene ID          | 400                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 62-181 of human ARL1 (NP_001168.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P40616                                                                                                           |
| Protein Weight   | 20kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |