

Anti-Arl13B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q3Sxy8

Item Number: CM0025893

Introduction

Rabbit polyclonal antibody against human ARL13B (JBTS8/ARL2L1), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying cilia, neural development, and Joubert syndrome. 49 kDa predicted weight.

[Learn More](#)

Parameter	Description
Product Name	ARL13B Rabbit pAb
Remarks/Alias	JBTS8; ARL2L1; ARL13B
Species	Human
GeneID	200894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ARL13B (NP_878899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3SXY8
Protein Weight	49kDa
Shipping	Ice bag