

Anti-Arl2 Polyclonal Antibody For Wb, Elisa - P36404

Item Number: CM0029853



Introduction

Anti-ARL2 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacting with Mouse and Rat. Recognizes human ARL2 (P36404). Ideal for cell cycle, mitochondrial, and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	ARL2 Rabbit pAb
Remarks/Alias	ARFL2; MRCS1; ARL2
Species	Human
GeneID (Human)	402
GeneID	402
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human ARL2 (NP_001658.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P36404
Protein Weight	21kDa
Shipping	Ice bag