

Anti-Arl3 Polyclonal Antibody For Wb, If/Icc, Elisa - P36405

Item Number: CM0022794

Introduction

High-quality rabbit polyclonal ARL3 antibody for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Recombinant immunogen, 20kDa target. Ideal for cilia, cell cycle, and protein trafficking research.

[Learn More](#)

Parameter	Value
Product Name	ARL3 Rabbit pAb
Remarks/Alias	RP83; ARFL3; JBTS35; ARL3
Species	Human
GeneID (Human)	403
GeneID	403
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human ARL3 (NP_004302.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36405
Protein Weight	20kDa
Shipping	Ice bag