

Anti-Arl4A Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P40617

Item Number: CM0019054

Introduction

Rabbit polyclonal antibody against human ARL4A, validated for IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat ARL4A protein. Suitable for studying small GTPase-mediated membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	ARL4A Rabbit pAb
Remarks/Alias	ARL4; ARL4A
Species	Human
Gene ID (Human)	10124
Gene ID	10124
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 121-200 of human ARL4A (NP_001032241.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40617
Protein Weight	23kDa
Shipping	Ice bag