

Anti-Arl4D Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49703

Item Number: CM0024094

Introduction

High-quality rabbit polyclonal antibody targeting ARL4D. Suitable for WB, IHC-P, and ELISA applications with cross-reactivity to human, mouse, and rat. Supports research into GTPase-mediated cellular processes.

[Learn More](#)

Parameter	Value
Product Name	ARL4D Rabbit pAb
Remarks/Alias	ARF4L; ARL4D
Species	Human
GenelD (Human)	379
GenelD	379
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 92-201 of human ARL4D (NP_001652.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49703
Protein Weight	22kDa
Shipping	Ice bag