


Anti-Arl6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H0F7

Item Number: CM0026039



Introduction

Rabbit polyclonal antibody targeting human ARL6 (BBS3, NP_115522.1). Validated for WB, IHC-P, ELISA; cross-reacts with mouse. Ideal for ciliary trafficking and Bardet-Biedl syndrome studies. ~21 kDa.

[Learn More](#)

Parameter	Value
Product Name	ARL6 Rabbit pAb
Remarks/Alias	BBS3; RP55; ARL6
Species	Human
Gene ID (Human)	84100
Gene ID	84100
Immunogen	Recombinant fusion protein containing amino acids 1-186 of human ARL6 (NP_115522.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0F7
Protein Weight	21kDa
Shipping	Ice bag