

Anti-Arl6Ip1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15041

Item Number: CM0023584

Introduction

Rabbit polyclonal antibody against ARL6IP1, validated for Western blot, IHC-P, and ELISA. Suitable for human, mouse, and rat samples. Useful for apoptosis, ER tubule, and glutamate transport research.

[Learn More](#)

Parameter	Specification
Product Name	ARL6IP1 Rabbit pAb
Remarks/Alias	AIP1; ARMER; SPG61; ARL6IP; ARL6IP1
Species	Human
GeneID (Human)	23204
GeneID	23204
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ARL6IP1 (NP_055976.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15041
Protein Weight	23kDa
Shipping	Ice bag