

Anti-Arl6Ip4 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q66Pj3

Item Number: CM0020951

Introduction

Rabbit polyclonal antibody targeting human ARL6IP4 (SRp25/SRp37), validated for IHC-P and ELISA. Ideal for alternative splicing studies and HSV-1 host-factor research. Cross-reacts with human. Protein weight 26kDa.

[Learn More](#)

Parameter	Value
Product Name	ARL6IP4 Rabbit pAb
Remarks/Alias	SR-25; SRp25; SFRS20; SRp37; ARL6IP4
Species	Human
GenelD (Human)	51329
GenelD	51329
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 313-402 of human ARL6IP4 (NP_057722.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q66Pj3
Protein Weight	26kDa
Shipping	Ice bag