

Anti-Arl6Ip5 Rabbit Polyclonal Antibody For Wb And Elisa - 075915

Item Number: CM0017612

Introduction

Rabbit polyclonal antibody against ARL6IP5 (PRAF3, JWA) for Western blot and ELISA applications. Detects human and mouse ARL6IP5, a regulator of glutamate and taurine transport. Uniprot O75915.

[Learn More](#)

Parameter	Value
Product Name	ARL6IP5 Rabbit pAb
Remarks/Alias	JWA; jmx; hp22; PRAF3; Yip6b; DERP11; HSPC127; addicisin; GTRAP3-18; ARL6IP5
Species	Human
Gene ID	10550
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human ARL6IP5 (NP_006398.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O75915
Protein Weight	22kDa
Shipping	Ice bag