


Anti-Arl8B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nvj2

Item Number: CM0025568



Introduction

High-quality rabbit polyclonal antibody specific to ARL8B. Suitable for WB, IF/ICC, and ELISA. Recognizes human and mouse ARL8B. Key tool for lysosomal trafficking, immune cytotoxicity, and antimicrobial pathway research.

[Learn More](#)

Parameter	Details
Product Name	ARL8B Rabbit pAb
Remarks/Alias	Gie1; ARL10C; ARL8B
Species	Human
Gene ID (Human)	55207
Gene ID	55207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-186 of human ARL8B (NP_060654.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NVJ2
Protein Weight	22kDa
Shipping	Ice bag