

Anti-Armc10 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N2F6

Item Number: CM0023943

Introduction

High-quality rabbit polyclonal antibody targeting human ARMC10, validated for Western blot and ELISA. Ideal for studying cell survival, growth regulation, and p53 signaling in various tissues.

[Learn More](#)

Parameter	Value
Product Name	ARMC10 Rabbit pAb
Remarks/Alias	SVH, PNAS112, PNAS-112, PSEC0198, ARMC10
Species	Human
GeneID (human)	83787
GeneID	83787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 169-308 of human ARMC10 (NP_001154481.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N2F6
Protein Weight	38kDa
Shipping	Ice bag